

BETA CELL KIT

<u>Item #</u>	<u>Item Name</u>	<u>List Price</u>
<u>0006-499</u>	Beta Cell Kit with crystal/sample retainer, cell holder and assorted spacers, without windows	\$123.00
<u>0006-498</u>	Beta Cell only crystal/sample retainer and assorted spacers	\$98.10
<u>0006-500</u>	Beta Cell Holder	\$55.10

The Beta Cell Kit is a convenient, low cost, easy to use demountable cell. This cell is used for qualitative analysis of liquids and Nujol mulls. Assembly is accomplished by placing the sample between two crystal windows separated by a spacer into the hexagonal holder and turning the single hexagonal threaded retainer until the cell is sealed.

The Beta Cell Kit uses standard 25mmx4mm crystal windows and spacers. The aperture of the cell is 19mm. The cell comes complete with a back plate, crystal/sample retainer and assorted spacers, but without window

INTERNATIONAL CRYSTAL LABORATORIES

11 Erie Street, Garfield, NJ 07026

www.internationalcrystal.net

iclmail@internationalcrystal.net

Phone: 973-478-8944.

